

Cell IDx 'RightON' Fluor and Biotin Antibody Labeling Kits

Description: Ultra-rapid 'mix and use' kits to incorporate optimal range of fluor or biotin on 50-100 µg of antibody. Optional final purification step included.

Cat. No.:

Labels	Cat#	Absorbance (nm)	Emission (nm)	Extinction Coefficient (M ⁻¹ cm ⁻¹)
Biotin-Peg4	ROA-XXX			
CL490	ROC-XXX	491	515	73,000
CL550	ROD-XXX	550	575	150,000
CL594	ROE-XXX	594	615	92,000
CL650	ROF-XXX	655	676	250,000
FITC (FAM)	ROK-XXX	488	490	70,000

Application: The Cell IDx 'RightON' fluor and biotin antibody labeling kits have been developed to modify 50-100 µg of purified antibody @ 0.5-1.0 mg/mL incorporating the optimal range of fluor or biotin in 10 min. Ultra-rapid, 'mix and use' protocol.

- **The kit also includes a ready-to-use spin filter to rapidly remove >80% of excess dye or biotin from final conjugate, with minimal loss, if desired.**
- All components are included in each kit
- Available labels include: biotin, FITC (FAM), CL405, CL490, CL550, CL594, CL650
- The only equipment needed are a pipette and a microcentrifuge for the optional excess dye removal step if desired
- Please ensure your antibodies for labeling are at 0.5-1.1mg/ml in buffer free of extraneous proteins such as BSA, serum or gelatin. Antibody formulations containing azide, trehalose or other sugars are acceptable. No tris buffer.

Components:

Kit Part	Quantity	Description
A	1	Modifier Tube
B	1	CellLight labeling reagent with lyophilized dye or Peg4-biotin for modification of 50-100 µg protein at 0.5-1.0 mg/mL. NOTE: IMMEDIATELY place and store biotin/dye @<-20°C
C	1	Red-capped desalting column to buffer exchange fluor- or biotin-modified antibody product into PBS + 15 mM sodium azide with collection tubes NOTE: Store columns @4°C
D	1	Final Buffer Tube

Storage: Fluor and biotin tubes @ -20°C, other components @ 2-8°C
Shelf-Life: 12 months

Protocol:

- Optimized for 50-100 µg of antibody @ 0.5-1.0 mg/mL.
- **Antibody needs to be protein carrier free, i.e., no gelatin, BSA, serum, etc.**
- **No Tris buffer**



Antibody Preparation

- 1) Add 1/10 volume of 10X modifier directly to your antibody tube (*i.e.* for 100 μ L of antibody, add 10 μ L of modifier).

Labeling with fluor or biotin

- 2) Add antibody solution to the biotin/dye tube, mix thoroughly with pipette or vortex. Incubate in the dark at room temperature for 10 minutes.

Final buffer step

- 3) After the 10-minute incubation is complete, add 90 μ L final buffer diluent to biotin/fluor-labeled antibody, vortex, and it is now ready to use in your assay!

Final antibody concentration is ~0.25-0.5 mg/mL, depending on starting concentration.

Notes

- For **biotin conjugation** it is highly recommended for optimal performance to proceed to the PURIFICATION STEP below following the 10-minute incubation with final buffer.

Optional Purification Step

- (a) PBS/azide desalting column preparation - Break off bottom of RED-capped column and loosen the RED cap, place column in a collection tube, then place in microcentrifuge and spin for 2 minutes @1500 g. Discard flow through.
- (b) Transfer probe labeled mixture to the RED-capped column, place column in a new collection tube and centrifuge for 2 minutes @1500 g.
- (c) Transfer desalted fluor labeled antibody to a new capped tube for 4^oC storage and you now have a dye/biotin free, labeled sample ready to use in your assay!

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*