

Absorption & Emission of CellLight Fluors Conjugated to Anti-hapten Detector Antibodies

Catalog #	Hapten	Fluor	Absorbance (nm)	Emission (nm)	Extinction Coefficient (M ⁻¹ cm ⁻¹)
UTK-490	CH015	CL490	491	515	73,000
UTK-550	CH016	CL550	550	575	150,000
UTK-595	CH019	CL595	594	615	92,000
UTK-650	CH021	CL650	655	676	250,000
UTK-750	CH019	CL750	759	780	240,000