



Product Specifications

UltraTag HRP1 detection antibody - UTCX-HRP1

Description:	Supplementary HRP labeled anti-UT016 detector antibody, 15 ug to stain 15-20 slides by IHC
Cat. No.:	UTCX-HRP1
Application:	Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC). Supplementary HRP labeled anti-UT016 detector antibody for use with UltraTag kits.
Components:	HRP labeled anti-UT016 detector antibody. 30ul at 0.5mg/ml. Suggested dilution 1/100.
Storage:	2-8°C
Shelf-Life:	Not determined.
Protocol:	LINK
Notes:	

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*