



---

## Product Specifications

### ***UltraTag AP1 detection antibody - UTCX-AP1***

Description:	Supplementary AP labeled anti-UT015 detector antibody, 15 ug to stain 15-20 slides by IHC
Cat. No.:	UTCX-AP1
Application:	Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC). Supplementary AP labeled anti-UT015 detector antibody for use with UltraTag kits
Components:	AP labeled anti-UT015 detector antibody. 30ul at 0.5mg/ml. Suggested dilution 1/100.
Storage:	2-8°C
Shelf-Life:	Not determined.
Protocol:	
Notes:	

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*