




6197 Cornerstone Court E, Ste 102
San Diego, CA 92121

Product Specifications

AlkPhos-anti-CH015 Conjugate

Description: Alkaline phosphatase (AlkPhos)-anti-CH015 antibody conjugate 15 µg to stain 15-20 slides, cells or spheroids.
Cat. No.: YEK-001
Application: Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC)
Components: AlkPhos-anti-CH015 antibody conjugate 15 µg
Storage: 2-8°C
Shelf-Life: Not determined.
Protocol: [LINK](#)

Notes:



(A) CD163-CH015 (2.5 µg/mL), anti-CH0163-AlkPhos (5.0 µg/mL) (B) UltraTag multiplex **Ki67/CD4/CD163**; clones Ki67 (AE1/AE3) CD163 (EP324), CD4 (EP204)

Disclaimer: *For In vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*

Version 1.0: Revision date 2.24.20