



Product Specifications

UltraTag HRP2 detection antibody - UTCX-HRP2

Description:	Supplementary HRP labeled anti-UT021 detector antibody, 15 ug to stain 15-20 slides by IHC
Cat. No.:	UTCX-HRP2
Application:	Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC). Supplementary HRP labeled anti-UT021 detector antibody for use with UltraTag kits
Components:	HRP labeled anti-UT021 detector antibody. 30ul at 0.5mg/ml. Suggested dilution 1/100.
Storage:	2-8°C
Shelf-Life:	12 months
Protocol:	
Notes:	

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*