



6197 Cornerstone Court, Ste 102
San Diego, CA 92121

Product Specifications

UltraPolymer Goat anti-Rabbit IgG (H&L) – HRP

| | |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description: | UltraPolymer Goat anti-Rabbit IgG (H&L) conjugated to HRP, affinity purified, min x w/ human, mouse or rat serum |
| Cat. No.: | 2RH-1000 |
| Amount: | 1000 mL |
| Concentration: | Ready to Use |
| Conjugate: | HRP |
| Applications: | IHC-P – optimized for use with autostainer or manual staining |
| Form: | Liquid |
| Buffer: | TBS, 1% (w/v) BSA, 0.1% Proclin 150, 10% heat-inactivated animal serum |
| Purification: | Affinity purified using solid phase Rabbit IgG |
| Purity: | Affinity purified antibody is > 95% based on SDS-PAGE |
| Host: | Goat |
| Clonality: | Polyclonal |
| Immunogen: | Purified Rabbit IgG, whole molecule |
| Storage: | 2-8 °C |
| Shelf-Life: | 1 year from received date when stored at 2-8 °C. |
| Specificity: | Based on IEP, this immunopolymer antibody reacts with heavy (γ) chains on rabbit IgG, light chains on all rabbit immunoglobulins |
| Cross-Reactivity: | Based on IEP, no cross-reactivity is observed to non-immunoglobulin rabbit serum proteins, serum proteins from human, mouse or rat |
| Country of Origin: | Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian. |
| Disclaimer: | <i>For In vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, inc. Liability of Cell IDx is limited to the purchase price of the product.</i> |

Tonsil stained with rabbit anti-human CD4 (clone EPR6855, 0.4ug/ml) and detected with UltraPolymer goat anti-mouse-HRP (RTU).
Staining on the BOND RX autostainer
Antigen retrieval: ER2, 20 min
Primary antibody: 30 min incubation
Secondary antibody: 15 min incubation

