



6197 Cornerstone Court E, Ste 102
San Diego, CA 92121

Product Specifications

UltraPolymer Donkey anti-Sheep IgG (H&L) – HRP

Description:	UltraPolymer Donkey anti-Sheep IgG (H&L) conjugated to HRP, affinity purified, min x w/ human, mouse or rabbit IgG
Cat. No.:	SH-050
Amount:	50ml
Conjugate:	HRP
Form:	Liquid
Buffer:	PBS, 1% BSA & 0.1%proclin 150
Purification:	Affinity purified using solid phase Sheep IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Donkey
Clonality:	Polyclonal
Immunogen:	Purified Sheep IgG, whole molecule
Storage:	2-8 °C
Shelf-Life:	1 year
Specificity:	Based on IEP, this immunopolymer antibody reacts with heavy (γ) chains on sheep IgG, light chains on all sheep immunoglobulins
Cross-Reactivity:	Based on IEP, no cross-reactivity is observed to non-immunoglobulin sheep serum proteins, IgG from human, mouse or rabbit
Country of Origin:	Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Disclaimer:	<i>For In vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, inc. Liability of Cell IDx is limited to the purchase price of the product.</i>