



6197 Cornerstone Court E, Ste 102  
San Diego, CA 92121

---

## Product Specifications

### UltraPolymer Donkey anti-Sheep IgG (H&L) – HRP

|                    |   |
|--------------------|---|
| Description:       | UltraPolymer Donkey anti-Sheep IgG (H&L) conjugated to HRP, affinity purified, min x w/ human, mouse or rabbit IgG  |
| Cat. No.:          | SH-100  |
| Amount:            | 1000ml  |
| Conjugate:         | HRP   |
| Form:              | Liquid  |
| Buffer:            | PBS, 1% BSA & 0.1%proclin 150   |
| Purification:      | Affinity purified using solid phase Sheep IgG   |
| Purity:            | Affinity purified antibody is > 95% based on SDS-PAGE   |
| Host:              | Donkey  |
| Clonality:         | Polyclonal  |
| Immunogen:         | Purified Sheep IgG, whole molecule  |
| Storage:           | 2-8 °C  |
| Shelf-Life:        | 1 year  |
| Specificity:       | Based on IEP, this immunopolymer antibody reacts with heavy ( $\gamma$ ) chains on sheep IgG, light chains on all sheep immunoglobulins   |
| Cross-Reactivity:  | Based on IEP, no cross-reactivity is observed to non-immunoglobulin sheep serum proteins, IgG from human, mouse or rabbit   |
| Country of Origin: | Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.  |
| Disclaimer:        | <i>For In vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, inc. Liability of Cell IDx is limited to the purchase price of the product.</i> |