



6197 Cornerstone Court E, Ste 102
San Diego, CA 92121

Product Specifications

UltraPolymer Goat anti-Mouse IgG (H&L) – HRP

Description:	UltraPolymer Goat anti-Mouse IgG (H&L) conjugated to HRP, affinity purified, min x w/ bovine, horse, human, pig or rabbit serum proteins
Cat. No.:	MH-100
Amount:	1000ml prediluted
Conjugate:	HRP
Form:	Liquid
Buffer:	PBS, 1% BSA & 0.1%proclin 150
Purification:	Affinity purified using solid phase Mouse IgG
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Purified Mouse IgG, whole molecule
Storage:	2-8 °C
Shelf-Life:	1 year
Specificity:	Based on IEP, this immunopolymer antibody reacts with heavy (γ) chains on mouse IgG, light chains on all mouse immunoglobulins
Cross-Reactivity:	Based on IEP, no cross-reactivity is observed to non-immunoglobulin mouse serum proteins, serum proteins from bovine, horse, human, pig or rabbit
Country of Origin:	Serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Disclaimer:	<i>For In vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, inc. Liability of Cell IDx is limited to the purchase price of the product.</i>

Triple Positive Breast Cancer Tissue
stained with Ki-67 and detected with
UltraPolymer goat anti-mouse HRP

