



Product Specifications

UltraTag 550 Fluorescent Detection Antibody - UTFX-550

Description:	Supplementary CL550 labeled anti-UT016 Tag detector antibody, 15 ug to stain 15-20 slides by IHC
Cat. No.:	UTFX-550
Application:	Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC). Supplementary CL550 labeled anti-UT016 detector antibody for use with UltraTag kits
Components:	CL550 labeled anti-UT016 detector antibody. 30ul at 0.5mg/ml. Suggested dilution 1/100.
Storage:	2-8°C
Shelf-Life:	Not determined.
Protocol:	
Notes:	

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*