



6197 Cornerstone Court E, Ste 102
San Diego, CA 92121

Product Specifications

CL550-anti-CH014

Description: CL550-anti-CH014 antibody 15 µg to stain 15-20 slides, cells or spheroids.
Cat. No.: UUH-550
Application: Multiplex immunoassays including fluorescent immunohistochemistry (IHC), immunocytochemistry (ICC) and ELISpot


α -Tag Probe	Fluor	Absorbance (nm)	Emission (nm)	Extinction Coefficient (M ⁻¹ cm ⁻¹)	Recommended Exposure Time (ms)
CH-014	CL-550	550	575	150,000	100-300

Suitable for spectral imaging.

Components: CL550-anti-CH014 antibody 15 µg

Storage: 2-8°C
Shelf-Life: Not determined.
Protocol: [LINK](#)

Notes:



(A) CD4-CH014/anti-CH014-CL550, (B) UltraTag multiplex CD163/CD4/CD68/LAG3; primary and anti-hapten antibody staining @ 5 µg/mL; clones CD163 (EPR19518), CD4 (EP204), CD68 (68/684)/LAG3 (294)

Disclaimer: *For In vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*

Version 1.0: Revision date 2.17.20