



6197 Cornerstone Court E, Ste 102
San Diego, CA 92121

Product Specifications

CL595-anti-CH019

Description: CL595-anti-CH019 antibody 15 µg to stain 15-20 slides, cells or spheroids.
Cat. No.: UTK-595
Application: Multiplex immunoassays including fluorescent immunohistochemistry (IHC), immunocytochemistry (ICC) and ELISpot

α -Tag Probe	Fluor	Absorbance (nm)	Emission (nm)	Extinction Coefficient (M ⁻¹ cm ⁻¹)	Recommended Exposure Time (ms)
CH019	CL595	594	615	92,000	400-800

Suitable for spectral imaging.

Components: CL595-anti-CH019 antibody 15 µg

Storage: 2-8°C
Shelf-Life: Not determined.
Protocol: [LINK](#)

Notes:

Disclaimer: *For In vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*

Version 1.0: Revision date 2.17.20