



---

## Product Specifications

### ***UltraTag 650 Fluorescent Detection Antibody - UTFX-650***

Description:	Supplementary CL650 labeled anti-UT021 Tag detector antibody, 15 ug to stain 15-20 slides by IHC
Cat. No.:	UTFX-650
Application:	Multiplex chromogenic immunoassays including chromogenic immunohistochemistry (IHC). Supplementary CL650 labeled anti-UT021 detector antibody for use with UltraTag kits
Components:	CL650 labeled anti-UT021 detector antibody. 30ul at 0.5mg/ml. Suggested dilution 1/100.
Storage:	2-8°C
Shelf-Life:	Not determined.
Protocol:	
Notes:	

Disclaimer: *For in vitro Research Use Only. Not for diagnostic or therapeutic use. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Cell IDx, Inc. Product may not be resold or modified for resale without prior written approval of Cell IDx, Inc.*